

CERTIFICATE OF ANALYSIS**Plasma/Serum Products****FT3 Serum**

Source Material:	Serum
Expiry Date:	None assigned. Users should determine suitability in their own system or refer to any known collection date.
<hr/>	
Product Code:	H625-2
Antigen/Marker present:	FT3
Species:	Human
Format:	Liquid
Storage Temperature:	≤-15°C. Minimise freeze-thaw cycles
Shipping Temperature:	≤-15°C recommended
Intended Use:	For performance analysis and quality assurance, research and further manufacturing

Health & Safety:

Infectious Disease Tests Performed:	The material has not been tested for the presence of infectious agents (other than under "Test Results Data"). See under Precaution for guidance.
Chemicals:	For further information please request a Safety Data Sheet.
Precaution:	Potentially hazardous. No test can guarantee absence of infective agents. Handle as hazardous using appropriate Personal Protective Equipment. Refer to Safety Data Sheet (available on request) for safety information and precautions.

Test Results Data:

Test performed: Labcorp FT3 Test #: 010389 (ng/ml)

Test result: See under "Lot/Donor information"


Lot/Donor information:

Lot #	Date of Collection	Donor Age	Donor Gender	Labcorp FT3 Test #: 010389 (ng/ml)
SE16212-A	06-Dec-2017	27	F	5.19
SE16213-A	18-Dec-2016	54	F	12.12
SE16214-A	17-Jul-2022	38	M	7.67
SE16215-A	07-Aug-2016	60	F	8.12
SE16216-A	06-Apr-2022	50	M	13.89
SE16217-A	11-Nov-2021	51	F	20.21
SE16218-A	14-Jul-2020	56	F	31.02
SE16219-A	30-Aug-2021	30	F	21.9
SE16220-A	23-Jan-2021	73	M	24.42
SE16221-A	28-Mar-2017	75	F	29.97

All Lots are remnants, from individual donors and not spiked or heat treated

In Logical Biological's lot numbering system, the lot number is shown before the dash (-). The letter after the dash is the "Tube Letter". Multiple tubes may be from exactly the same lot, e.g., PL001-A and PL001-B would be two tubes of exactly the same lot

Quality Assurance Approval:

Signature:	
Date:	13 th March 2023